

Starkey's Revolutionary Product Line

Like nothing else. From no one else.



Hearing Aid Styles

livioEdge^{AI}
livio^{AI}
livio



Feature	BTE R	ITC/HS R and ITE R	RIC R	RIC 312	micro RIC 312	BTE 13
Battery Type	Rechargeable Lithium-ion	Rechargeable Lithium-ion	Rechargeable Lithium-ion	312	312	13
Wireless Protocols	2.4 GHz	2.4 GHz	2.4 GHZ + NFMI	2.4 GHz + NFMI	2.4 GHz	2.4 GHZ + NFMI
Telecoil	•	• (ITE R only)	•			•
CROS System			•	•		•
Tinnitus Technology	•	•	•	•	•	•

* Livio 1000 only available in BTE R and BTE 13 styles

2.4 GHz ACCESSORIES



Remote Microphone +
Use as a remote microphone, or easily stream binaurally from four different audio sources: Bluetooth®, Telecoil, FM, or Line-In.



Remote
Our Remote control includes updated features to control memory and volume, mute hearing aids and turn other special features on and off.



Mini Remote Microphone
Enjoy one-on-one conversations in noisy environments with our small, easy-to-use Mini Remote Microphone. Clip it onto the clothing of the person being spoken to, or use it as a TV streamer by placing it near the speaker.

Standard Charger & Custom Charger

Both chargers hold enough power to provide portable charging three times. LED indicators let patients know when their hearing aids are fully charged, and hearing aids turn on automatically when removed from the charger.




Mini Turbo Charger

The Mini Turbo Charger is a pocket-size lithium-ion charging option providing a 3.5 hour charge in just 7 minutes or a full charge in under 3 hours. This portable charger will recharge a pair of RIC R or BTE R hearing aids up to four times without a power connection.



Hearing Performance & Healthable Technology Features

AI Features by Technology Tier



Edge Mode

Instantly conducts an AI-based analysis of the listening environments and makes smart, immediate adjustments



IntelliVoice

Thanks to the power of AI, this first-of-its-kind speech enhancement feature uses the wearer's iOS smartphone to deliver an unprecedented experience to those with a moderate loss of 50 dB HL or greater



Thrive Assistant

Patients use voice commands to get answers to questions, set reminders and control their hearing aids



Find My Phone

Using their hearing aids and Thrive Assistant, wearers are able to locate their lost or misplaced phone by having it ring



Fall Alert

Detects when a hearing aid wearer falls and sends an alert message to selected contacts



Translate

Patients instantly hear speech translated in their language right into their hearing aids



Transcribe

Conversations are transcribed into text that can be saved or shared



Reminders

Patients get in-ear and on-screen notifications for tasks like taking medications












Auto On/Off

Automatically and conveniently saves battery power when the patients hearing aids are not in use


*Reminders available for Livio AI 1200 as a Thrive In-App purchase.

The app may look slightly different depending on your phone and/or hearing aid model.

The Starkey logo, Starkey, Livio, Livio logo, Thrive, Hearing Reality and Hearing Care Anywhere are registered trademarks of Starkey Laboratories, Inc. The Thrive logo, Thrive Care, IntelliVoice and Healthable are trademarks of Starkey Laboratories, Inc. OrCam is a registered trademark of OrCam Technologies LTD. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks is under license. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

- 
Android and iOS Streaming
 Stream phone calls and media directly from not only Apple/iOS devices, but also selected Android phones with the Android 10 operating system
- 
Hearing Care Anywhere
 Allows Livio Edge AI, Livio AI or Livio patients to receive remote adjustments to their hearing aids without having to stop into the office
- 
Tap Control
 Enables you to easily stop or start audio streaming, activate Edge Mode or access the Thrive Assistant by simply tapping your ear twice
- 
Self Check
 Provides a quick, convenient way for patients to analyse their hearing aid system's performance
- 
Body Tracking
 Tracks daily steps, measures movement and monitors more vigorous physical activity
- 
Brain Tracking
 Monitors the brain-health benefits of wearing hearing aids by tracking things like social engagement
- 
Find My Hearing Aids
 Helps patients easily locate misplaced hearing aids providing a location and a timestamp
- 
Personalised Memories
 Create multiple custom memories including geotagged memories like, 'coffee shop' or 'home'
- 
Remote Control
 Patients can change volume and switch memories easily using a tool in the Thrive app.

Two powerful apps that work together















































































































thriveCare 
 For the people who care about them

Thrive Care is exclusive to Livio Edge AI.

Features Overview

HEARING REALITY PRO

	Livio Edge AI Premium 2400	Livio Premium 2400	Livio AI/Livio Advanced 2000	Livio AI/Livio Select 1600	Livio AI/Livio Basic 1200	Livio Economy 1000
Sound Imaging: Channels Bands	24	24	20	16	12	10
Speech Optimization						
Music Optimization:						
Music Adaptation						
E2E Music Adaptation						
Ear-to-Ear Technology:						
E2E Wind Noise Management						
E2E Machine Noise Adaptation						
E2E Directionality						
E2E Phone Streaming						
Edge Mode						
IntelliVoice						
Sound Manager:						
Auto Music						
Quiet						
Spatial Speech Enhancement						
Machine Noise	Up to 20 dB of reduction	Up to 20 dB of reduction	Up to 10 dB of reduction	Up to 5 dB of reduction	Up to 5 dB of reduction	Up to 5 dB of reduction
Speech in Noise	Up to 20 dB of reduction	Up to 20 dB of reduction	Up to 8 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction
Transient Noise Reduction	Up to 15 dB of reduction	Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction
Wind	Up to 30 dB of reduction	Up to 30 dB of reduction	Up to 15 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction
Directional Processing:						
Immersion						
Adaptive						
Dynamic						
Directional						
Feedback Management						
Frequency Lowering						
Hearing Care Anywhere						
Self Check						
REM Target Match						
2.4 GHz Accessory Compatibility						

The app may look slightly different depending on your phone and/or hearing aid model.

The Starkey logo, Starkey, Livio, Livio logo, Thrive, Hearing Reality and Hearing Care Anywhere are registered trademarks of Starkey Laboratories, Inc. The Thrive logo, Thrive Care, IntelliVoice and Healthable are trademarks of Starkey Laboratories, Inc.

©2021 Starkey. All Rights Reserved. FLYR_202107_AU