



Edge Mode+

Giving patients control of what they want to hear when it matters most

How it helps

Edge Mode+ is Starkey’s second generation intelligent, on-demand feature that enables instant, powerful speech audibility and noise management adjustments when patients need it most. Edge Mode+ gives patients the ability to automatically select when they need more help, eliminating the need for manual programs.

How it works

Once activated, Edge Mode+ instantly conducts an acoustic analysis of the listener’s environment and harnesses Deep Neural Network computing to make automatic adjustments to gain, noise management, and directionality. With Genesis AI, patients can now select their listening intent—enhance speech or reduce noise—to further personalize their Edge Mode+ setting. When the patient moves out of that environment, they simply deselect Edge Mode+ and enjoy the automatic noise management adjustments that the hearing aids provide.




Figure 1: Edge Mode+ control options in the My Starkey mobile app

Perceived Benefit





Figure 2: Clinical testing responses that indicated a benefit of Edge Mode+ (n=44) and selecting listening intent (n=35).

Where to find it



Patient Feedback



Edge Mode+ is impressive. It worked so good. It was almost unsettling at first because I felt like I wasn’t hearing well, and then what I realized is I wasn’t hearing the noise well, but I was understanding voices perfectly.

