

GENESIS^{AI}



PRODUCT GUIDE

2024



GENESIS AI PRODUCT GUIDE 2024

	Performance				
Premium 24	mRIC R	RIC RT	RIC 312	ITC R	ITE R
Advanced 20	mRIC R	RIC RT	RIC 312	ITC R	ITE R
Select 16	mRIC R	RIC RT	RIC 312	ITC R	ITE R

Discreet Design			
IIC NW	CIC NW	CIC	mRIC R
IIC NW	CIC NW	CIC	mRIC R
IIC NW	CIC NW	CIC	mRIC R

Power Custom Cased SnapFit 2.0 Receiver		
mRIC R	RIC RT	RIC 312
mRIC R	RIC RT	RIC 312
mRIC R	RIC RT	RIC 312

CROS	
RIC RT	RIC 312
RIC RT	RIC 312
RIC RT	RIC 312

Available Colors

- | | | | | | |
|-----------------|---------------|------------------|--------------|--------------|-------------|
| Standard Colors | White | Faceplate Colors | Pink | Shell Colors | Clear |
| | Tech Black | | Light Brown | | Light Brown |
| | Graphite Gray | | Medium Brown | | Pink |
| | Silver | | Chestnut | | Blue/Red |
| | Chestnut | | Dark Brown | | Black |
| | Caramel | | Black | | |
| | Beige | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

StarLink Wireless Accessories

<p>StarLink Table Microphone Make group settings easy with eight built-in microphones that stream the primary speaker's voice directly to your patient's hearing aids.</p>	<p>StarLink Remote Microphone + Allows patients to easily stream from a variety of audio sources. This simple-to-use remote microphone + offers the most capabilities available in any accessory.</p>	<p>StarLink Mini Remote Microphone Enables easy audio streaming and one-on-one conversations in noisy environments.</p>
<p>StarLink Remote Control 2.0 Easily control memory and volume, mute hearing aids and turn other special features on and off.</p>	<p>StarLink TV Streamer Allows patients to easily stream audio from their TV or other electronic audio source directly to their hearing aids.</p>	<p>StarLink Charger 2.0 All chargers come with on-board turbo charging. Turbo charging automatically occurs when the hearing aid battery is low. The hearing aid will rapidly charge, providing approximately 3 hours of use time with only 10 minutes of charging. Premium Charger with on-board battery available for on-the-go charging.</p>

Genesis AI introduces an all-new design for maximum quality and durability. Every rechargeable Genesis AI product is waterproof* thanks to our Pro8 HydraShield™ technology.

Waterproof

Humidity Test	Salt Mist Test	Rain Spray Test	Targeted Sweat Drip Test**	Extended Reliability Testing**	Salt Brine Immersion Test**
IP68+					

*Up to 1m. **RIC Genesis AI products only.

Genesis AI Key Features

Our latest technology provides superior sound quality plus advanced features and capabilities.

	Premium 24	Advanced 20	Select 16
Voice AI	●		
Translate	●		
Personal Activity Monitoring	●		
Edge Mode + Automatic Adaptation	●		
Edge Mode+ with user control	●		
Feedback Canceller – Motion-Based Optimisation	●		
Directional Microphone – Motion-Based Optimisation	●	●	
Multiflex Tinnitus Pro	●	●	
Transcribe	●	●	
Intelligent Reminders	●	●	
Edge Mode+ – Personalisation	●	●	
Edge Mode + MyStarkey app access	●	●	●
Smart Assistant with Voice	●	●	●
Find My Phone	●	●	●
Hear Share app compatibility	●	●	●
Engagement Tracking	●	●	●
Activity Tracking	●	●	●
Fall Alert	●	●	●
2-Way Audio	●	●	●
Basic Reminders	●	●	●
Self Check	●	●	●
Auto Sleep	●	●	●
Tap Control	●	●	●
Tinnitus with Tilt Control	●	●	●
TeleHear Remote Programming (Live Sessions & Asynchronous)	●	●	●
REM Target Match	●	●	●

App compatibility

All Genesis AI wireless styles are compatible with our My Starkey patient app and Hear Share caregiver app.



My Starkey patient app

Gives patients control over their hearing aids, plus the ability to get tips, track their health, and access intelligent features designed to simplify their lives.



Hear Share caregiver app

Provides peace of mind for caregivers and more independence for hearing aid wearers.

Pairing Genesis AI with Starkey's popular apps gives patients a long list of extra features, including:

Edge Mode+

Uses AI to optimize sound quality on-demand for either clearer speech or improved listening comfort.



Health monitoring

Uses the science of the ears to accurately count steps, track engagement and help wearers achieve well-being goals.



Fall detection

Starkey is still the only manufacturer with hearing aids that can detect falls and alert family and friends.



Starkey Neuro Processor

A processor with a mind of its own.

The processor is the brains of every hearing aid, and no processor is smaller or more powerful than our all-new Starkey Neuro Processor.

It's the industry's most advanced processor technology with a unique on-board Deep Neural Network (DNN) accelerator engine that mimics the cerebral cortex of the human brain.

Tech Tiers Channels Bands	Premium 24	Advanced 20	Select 16
Edge Mode+			
User Control	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Automatic Adaptation	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Personalization	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
App Control	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Speech Optimization			
Consonant Brightness	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Voice AI	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Sound Manager			
Ear-to-Ear Technology	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Quiet	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Machine Noise Reduction	<p>Up to 22 dB of reduction</p>	<p>Up to 12 dB of reduction</p>	<p>Up to 7 dB of reduction</p>
Speech in Noise Reduction	<p>Up to 22 dB of reduction</p>	<p>Up to 10 dB of reduction</p>	<p>Up to 8 dB of reduction</p>
Transient Noise Reduction	<p>Up to 15 dB of reduction</p>	<p>Up to 9 dB of reduction</p>	<p>Up to 6 dB of reduction</p>
Wind Reduction	<p>Up to 35 dB of reduction</p>	<p>Up to 19 dB of reduction</p>	<p>Up to 10 dB of reduction</p>
Directional Processing*			
Motion-Based Optimisation	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
Adaptive	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Dynamic	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>

*Default on all standard styles. Optional for ITC/HS R and ITE R styles.

