



Product Portfolio Guide

2025 | Australia


EDGE^{AI}

Tech Tier/Channels/Bands

24 / 20 / 16

AURACAST

CROS COMPATIBLE



- Waterproof*
- Wireless
- Rechargeable

*Up to 1m.




Available Colors




- Standard Colors**
- White
 - Tech Black
 - Graphite Gray
 - Silver
 - Chestnut
 - Caramel
 - Beige

- Custom Colors**
- | | |
|-------------------------|---------------------|
| Faceplate Colors | Shell Colors |
| Pink | Clear |
| Light Brown | Pink |
| Medium Brown | Light Brown |
| Chestnut | Blue/Red |
| Dark Brown | Black |
| Black | |

SnapFit 2.0 Stock and Custom Cased Receivers



Fitting Ranges



Visit StarkeyPro.com/EdgeAI for full technical data sheets.



Powered by the G2 Neuro Processor with the industry's only neural processing unit (NPU) integrated into the chip, and an all-new Bluetooth® Low Energy (LE) Audio streaming protocol, Starkey Edge AI's advanced processing is able to provide so much more than all-day battery life and waterproof protection – allowing patients to live easier, more connected lives.

For more Edge AI details, visit StarkeyPro.com.au























EDGE^{AI}

Tech Tiers | Channels | Bands

Premium
24

Advanced
20

Select
16

Edge Mode+			
Automatic	●	○	○
User Control	●	○	○
Personalization	●	●	○
App Control	●	●	●
Speech Optimization			
Consonant Brightness	●	●	●
Voice AI*	●	○	○
Sound Manager			
Quiet 	 Maximum expansion	 Moderate expansion	 Mild expansion
Machine Noise Reduction 	 Up to 22 dB of reduction	 Up to 12 dB of reduction	 Up to 7 dB of reduction
Speech in Noise Reduction 	 Up to 22 dB of reduction	 Up to 10 dB of reduction	 Up to 8 dB of reduction
Transient Noise Reduction 	 Up to 15 dB of reduction	 Up to 9 dB of reduction	 Up to 6 dB of reduction
Wind Reduction 	 Up to 35 dB of reduction	 Up to 19 dB of reduction	 Up to 10 dB of reduction
Additional Key Features			
Translate	●	○	○
Transcribe	●	●	○
Multiflex Tinnitus Pro with new tinnitus maskers	●	●	○
Smart Assistant & Find My Hearing Aids	●	●	●
2-Way Audio	●	●	●
Hearing Health & Activity Tracking	●	●	●
Fall Alert & Balance Assessment	●	●	●
Reminders & Self Check	●	●	●
Auto Sleep	●	●	●
Tap Control	●	●	●
TeleHear full remote fit	●	●	●
REM Target Match	●	●	●
Demo Stream	●	●	●
Waterproof up to 1m**	●	●	●
Apple Watch Controls*	●	●	●

Features noted above in blue are available within My Starkey app.

*For iOS users only. **Rechargeable styles only.

Wireless Accessories

ALL NEW



StarLink Edge TV Streamer
Part No.P00004316

AURACAST



StarLink Remote Control 2.0
Part No. P00002347

Charging Solutions

All chargers include
Turbo Charge



StarLink Premium Charger 2.0
w/ onboard battery

- mRIC R
- RIC RT

StarLink Charger 2.0***

- mRIC R
- RIC RT

***Default selection with rechargeable hearing aid order.



StarLink Premium Custom Charger 2.0
w/ onboard battery

- ITC/HS/ITE R


StarLink Custom Charger 2.0***

- ITC/HS/ITE R

Compatible Apps



My Starkey



Hear Share



Scan to check
phone compatibility



Tech Tier/Channels/Bands
24 / 20* / 16*



Waterproof**
Rechargeable

For more product details,
visit StarkeyPro.com.au



Faceplate Colors

- Pink
- Light Brown
- Medium Brown
- Chestnut*
- Dark Brown
- Black

Shell Colors

(IIC only available in Black)

- Clear
- Pink
- Light Brown
- Blue/Red
- Black

Charging Solutions

Both chargers include
Turbo Charge




StarLink Premium Custom Charger 2.0
w/ onboard battery

- CIC R NW

StarLink Custom Charger 2.0***

- CIC R NW

Fitting Ranges



Based on the industry's most advanced processor and powered by our industry-leading Neuro Sound Technology, Starkey Signature Series deliver the size, performance and personalized fit that patients can only get from you.

Tech Tiers Channels Bands	Premium 24	Advanced 20*	Select 16*
Edge Mode+ [CIC styles only]			
User Control	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Speech Optimization			
Consonant Brightness	<input checked="" type="radio"/> Maximum	<input checked="" type="radio"/> Mild	<input checked="" type="radio"/> Mild
Sound Manager			
Quiet	<input checked="" type="radio"/> Maximum expansion	<input checked="" type="radio"/> Moderate expansion	<input checked="" type="radio"/> Mild expansion
Machine Noise Reduction	<input checked="" type="radio"/> Up to 22 dB of reduction	<input checked="" type="radio"/> Up to 12 dB of reduction	<input checked="" type="radio"/> Up to 7 dB of reduction
Speech in Noise Reduction	<input checked="" type="radio"/> Up to 22 dB of reduction	<input checked="" type="radio"/> Up to 10 dB of reduction	<input checked="" type="radio"/> Up to 8 dB of reduction
Transient Noise Reduction	<input checked="" type="radio"/> Up to 15 dB of reduction	<input checked="" type="radio"/> Up to 9 dB of reduction	<input checked="" type="radio"/> Up to 6 dB of reduction
Wind Reduction	<input checked="" type="radio"/> Up to 35 dB of reduction	<input checked="" type="radio"/> Up to 19 dB of reduction	<input checked="" type="radio"/> Up to 10 dB of reduction
Additional Key Features			
Multiflex Tinnitus Pro with new tinnitus maskers	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
REM Target Match	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Waterproof**	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>

*CIC styles only. **Up to 1m. Rechargeable style only. ***Default selection with rechargeable hearing aid order.



Quality in... quality out

Your professional expertise is critical when selecting the best solution for a patient. Many factors and considerations come into play when determining if a patient is a custom hearing aid candidate, including:

- Audiogram
- Dexterity
- Ear anatomy and history
- Lifestyle

Ultimately, **the patient has come to you** to weigh all considerations—physical, audiological and patient priorities—and recommend a solution that delivers both success and satisfaction.

A **superior custom product fitting starts with a high quality ear impression**. Making a perfect impression requires proper technique, an attention to detail, and repetition.

As the custom hearing aid experts, we know a thing or two about making accurate, precision-tested impressions —and how good impressions can also be used to make earmolds, custom cases, hearing protection solutions and more.



Scan the QR code or go to StarkeyPro.com.au for support materials, how-to videos, and more.

Developed with help from Starkey's most experienced impression takers, these tools share best practice tips perfected from shooting thousands of impressions in our Center for Excellence.

GENESIS^{AI}

Tech Tier/Channels/Bands
24 / 20 / 16 / 12



Waterproof*
 Wireless
 Rechargeable
 *Up to 1m.



Available Colors




- Standard Colors**
- White
 - Tech Black
 - Graphite Gray
 - Silver
 - Chestnut
 - Caramel
 - Beige

- Custom Colors**
- Faceplate Colors**
- Pink
 - Light Brown
 - Medium Brown
 - Chestnut
 - Dark Brown
 - Black
- Shell Colors**
- Clear
 - Pink
 - Light Brown
 - Blue/Red
 - Black

SnapFit 2.0 Stock and Custom Cased Receivers



Fitting Ranges



Visit StarkeyPro.com for full technical data sheets.



Genesis AI marked a new beginning for hearing technology by introducing the industry to the most advanced processor technology with a unique onboard Deep Neural Network (DNN) accelerator. This Neuro Sound Technology in Genesis AI is designed to deliver sound so transparent and true to life, patients will forget they're wearing hearing aids.

For more product details, visit StarkeyPro.com.au



GENESIS^{AI}

Tech Tiers | Channels | Bands

Premium **24** Advanced **20** Select **16** Basic **12**

Edge Mode+		Premium 24	Advanced 20	Select 16	Basic 12
Automatic User Control Personalization	App Control	●	○	○	○
Speech Optimization					
Consonant Brightness	Voice AI*	●	○	○	○
Sound Manager					
Quiet		Maximum expansion	Moderate expansion	Mild expansion	Mild expansion
Machine Noise Reduction		Up to 22 dB of reduction	Up to 12 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction
Speech in Noise Reduction		Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction
Transient Noise Reduction		Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction
Wind Reduction		Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction	Up to 10 dB of reduction
Additional Key Features					
Translate	Transcribe	●	○	○	○
Multiflex Tinnitus Pro	Smart Assistant & Find My Hearing Aids	●	●	○	○
2-Way Audio	Hearing Health & Activity Tracking	●	●	●	●
Fall Alert	Reminders & Self Check	●	●	●	●
Auto Sleep	Tap Control	●	●	●	●
TeleHear full remote fit	REM Target Match	●	●	●	●
Waterproof up to 1m**	Apple Watch Controls*	●	●	●	●

Features noted above in blue are available within My Starkey app.

*For iOS users only. **Rechargeable styles only.

Wireless Accessories



StarLink Remote Control 2.0
Part No. P00002347



StarLink Mini Remote Microphone
Part No. 43121-000



StarLink Remote Microphone +
Part No. 43120-000



StarLink Table Microphone
Part No. 43123-610



StarLink TV Streamer
Part No. 43119-000

Charging Solutions

All chargers include
— Turbo Charge —



StarLink Premium Charger 2.0
w/ onboard battery

- mRIC R
- RIC RT

StarLink Charger 2.0***

- mRIC R
- RIC RT

***Default selection with rechargeable hearing aid order.



StarLink Premium Custom Charger 2.0
w/ onboard battery

- ITC/HS/ITE R


StarLink Custom Charger 2.0***

- ITC/HS/ITE R

Compatible Apps



My Starkey



HEAR Share



Scan to check
phone compatibility



evolv^{AI}

Tech Tiers
2400 / 2000 / 1600 / 1200 / 1000




- Wireless
- Rechargeable


Available Colors

- Champagne
- Black
- Slate
- Espresso
- Sterling
- Bright White with Sterling
- Bronze





BTE Stock Solutions



SnapFit Stock and Absolute Power Receivers



Fitting Ranges







































Visit StarkeyPro.com for full technical data sheets.

Evolv AI offers a complete line of BTE (Behind-The-Ear) styles loaded with innovative technology and clear sound that seamlessly connect your patients to the life they love.

For more product details, visit StarkeyPro.com.au



*RIC R not available in 1000 technology level

	Premium 2400 24	Advanced 2000 20	Select 1600 16	Basic 1200 12	Entry 1000 10
Edge Mode					
User Control	●	○	○	○	○
App Control	●	●	●	○	○
Speech Optimization					
Spatial Speech Enhancement	●	●	●	○	○
Ear-to-Ear Technology	●	●	●	○	○
Voice AI*	●	○	○	○	○
Sound Manager					
Auto Music 	 Strong	 Moderate	 Mild		
Quiet 	 Maximum expansion	 Moderate expansion	 Mild expansion	 Mild expansion	 Mild expansion
Machine Noise Reduction 	 Up to 22 dB of reduction	 Up to 12 dB of reduction	 Up to 7 dB of reduction	 Up to 7 dB of reduction	 Up to 7 dB of reduction
Speech in Noise Reduction 	 Up to 22 dB of reduction	 Up to 10 dB of reduction	 Up to 8 dB of reduction	 Up to 8 dB of reduction	 Up to 8 dB of reduction
Transient Noise Reduction 	 Up to 15 dB of reduction	 Up to 9 dB of reduction	 Up to 8 dB of reduction	 Up to 8 dB of reduction	 Up to 8 dB of reduction
Wind Reduction 	 Up to 32 dB of reduction	 Up to 17 dB of reduction	 Up to 9 dB of reduction	 Up to 9 dB of reduction	 Up to 9 dB of reduction
Additional Key Features					
Translate	●	○	○	○	○
Multiflex Tinnitus Pro	●	○	○	○	○
Transcribe	●	●	○	○	○
Thrive Assistant & Find My Hearing Aids	●	●	●	○	○
CROS System**	●	●	●	●	○
Engagement & Activity Tracking	●	●	●	●	○
Fall Alert	●	●	●	●	○
2-Way Audio	●	●	●	●	●
Reminders & Self Check	●	●	●	●	●
Auto On/Off	●	●	●	●	●
Tap Control	●	●	●	●	●
TeleHear full remote fit	●	●	●	●	●
REM Target Match	●	●	●	●	●

Features noted above in blue are available within Thrive app. *For iOS users only. ** BTE 13, RIC 312 & RIC R styles only.

Wireless Accessories



StarLink Remote Control 2.0
Part No. P00002347



StarLink Mini Remote Microphone
Part No. 43121-000



StarLink Remote Microphone +
Part No. 43120-000



StarLink Table Microphone
Part No. 43123-610



StarLink TV Streamer
Part No. 43119-000

Charging Solutions



Standard Premium Charger w/ onboard battery

- BTE R & RIC R

Standard Charger 2.0***


- BTE R & RIC R

***Default selection with rechargeable hearing aid order.

Compatible Apps



Thrive Hearing Control



Hear Share



Scan to check phone compatibility



Earmolds

As the leader in customization, Starkey understands offering your patients an expansive portfolio of hearing care solutions specific to their needs is key.

FIND ME:
eStore > Shop > Earmolds



Hear better. Live better.

Material Options

Digital SLS Acrylic	<ul style="list-style-type: none"> RIC BTE & Thin Tube Specialty 	<ul style="list-style-type: none"> One color 	Most Rigid 90 Shore
High Strength Silicone 60	<ul style="list-style-type: none"> RIC BTE & Thin Tube Specialty 	<ul style="list-style-type: none"> Swirl 3 colors Glitter option 	60 Shore
High Strength Silicone 40	<ul style="list-style-type: none"> RIC BTE & Thin Tube Specialty 	<ul style="list-style-type: none"> Swirl 3 colors Glitter option 	40 Shore
Floatable Silicone	<ul style="list-style-type: none"> RIC BTE & Thin Tube Specialty 	<ul style="list-style-type: none"> Swirl 3 colors 	25 Shore Softest

Venting Guide

This guide is applicable for all earmolds and earplugs. Vent size is limited by ear anatomy.

Vent Name	Approximate Size
Pressure Vent	0.8 mm
Small	1.2 mm
Medium	1.8 mm
Large	2.55 mm
Select-A-Vent	3.37 mm
Extra Large	As large as canal allows

Color Options

Our earmolds come in a variety of color options. **The latest color options can be found on the eStore.** Some colors appear differently in certain color families based on material type. Swirl combinations are also available using up to three colors.

Pro Tip: Adding White as one of the three colors within a color swirl combination will help separate and make the other color choices "pop."


Color combinations to avoid:

- 'Light' colors combined in Swirl or 'Light' colors combined with 'Clear'
- Similar colors combined in Swirl

RIC Solutions

RIC Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Available for both SnapFit 2.0 and SnapFit Receivers




Custom Cased or Absolute Power (AP) Receivers Digital SLS

Cable Lengths:
1, 2, 3, 4, 5

Gain:
L, M, P, UP






Receivers:
SnapFit 2.0
(ten-pin connector)

SnapFit
(six-pin connector)
Original
(two-pin connector)



*Only available in clear
Pro tip: When measuring the receiver wire for a Custom Cased or AP receiver, report the number that correlates with the top of the entrance to the ear canal.


Earplug Solutions (Hearing Conservation)



BTE Solutions

Thin Tube Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40


Tubing Lengths:











Traditional Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Tubing:

- 13 Regular
- 13 Heavy
- 13 Extra Heavy
- 13 Dry Tubing
- 13 Regular Snap Connect
- 13 Heavy Snap Connect




Specialty Solutions





Future *Forward*

By continuously pushing the edge of what's possible, Starkey will always be 10 years ahead—with groundbreaking innovations, award-winning products and services, and a **complete ecosystem of solutions** meticulously integrated to help you effortlessly provide better hearing for all.



Hearing is Our Concern.

Starkey logo, Starkey, TeleHear, Pro Fit, Starkey Edge, Edge Mode, Multiflex Tinnitus, StarLink, Starkey Signature Series, Neuro Sound Technology, Genesis, Hear Share, Evolv, Thrive, SoundGear logo, SoundGear, Phantom and Tunz are trademarks of Starkey Laboratories, Inc.

Hearing is Our Concern is a registered trademark of Starkey Laboratories, Inc.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by the Bluetooth SIG. Any use of such marks by Starkey is under license. Other trademarks and trade names are those of their respective owners.

Apple Watch is a trademark of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. "Made for iPod," "Made for iPhone" and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, iPad or Apple Watch may affect wireless performance. Google Play Store is a trademark of Google, LLC.